

Sante PACS Server

The Web Server Module

Copyright © 2026 [Santesoft LTD](#), all rights reserved

Table of Contents

Overview	2
Free Applications	2
Configure the HTTP Server	3
Enable SSL Encryption	4
Study Accessibility	4
How to use the server	5
Add Server User	6
Server Activity	9
Remote connection with an internet browser	12
Study retrieving by patients	17
Study retrieving with use of token	17
Token with fully qualified web address	19
Sante DICOM Web Viewer	20
The OHIF (Open Health Imaging Foundation) Viewer	21
The MIT License	21
The interface of the OHIF Viewer	21
Customize the web page appearance	22
Customize the logo image, the favicon, the web site's address and the popup alternative text of the logo image	23
Customize the header's text, text and background color	24
Customize the quote's text and background color	25
Customize the unread study's text and background color, and the subtitle's text and background color	26
Customize the read study's text and background color, and the subtitle's text and background color	27
Customize the button's text and background color, and the hover text and background color	28
Customize the footer's text, text and background color and the company's link	29

Overview

The Web Server module allows remote users to connect to the system and perform queries and retrieve studies stored in the DICOM Server database. The module uses the HTTP protocol for communication, and remote users must use a web browser such as **Microsoft Edge**, **Google Chrome**, or **Mozilla Firefox** to access the server and retrieve studies.


Free Applications

From the initial web page of the Web Server, the remote user can download and install two free applications designed for use with this server: **Sante Web Downloader** and **Sante DICOM Web Viewer**.



Sante Web Downloader allows the user to monitor the downloading process. When the download is complete, the application saves the received studies in a ZIP file. The user can configure **Sante Web Downloader** to automatically open the saved ZIP file with a preferred DICOM viewer, such as **Sante DICOM Viewer Pro** or **Sante DICOM Editor**. Note that the selected DICOM viewer must support opening ZIP files from the command line. For more information, please refer to the *Sante Web Downloader User Guide*. This application is particularly useful when a remote user prefers to view studies with a DICOM viewer other than **Sante DICOM Web Viewer**.

Sante DICOM Web Viewer also allows the user to monitor the downloading process, similar to **Sante Web Downloader**, but it additionally functions as a fully featured, free 2D/3D DICOM viewer. It provides a comprehensive set of tools, including measurement tools, MPR (Multi-Planar Reconstruction), and study/series comparison. The viewer is designed to operate smoothly in multi-modality and multi-manufacturer environments. For more information, please refer to the *Sante DICOM Web Viewer User Guide*.



Sante Web Downloader



Sante DICOM Web Viewer

Configure the HTTP Server

Select the menu command "**Network** ➔ **HTTP Server Setup**".

A dialog box appears, allowing the user to configure the HTTP Server module.

Internet IP Address (True IP)

The Internet IP address is automatically assigned by the Internet Service Provider (ISP). It is the address through which the HTTP Server communicates with the Internet. This address cannot be changed by the user.

Web IP Port

This is the port number used by the HTTP Server for Internet communication. If the Web Port value is set to **0**, the program does not open the port, and therefore the server cannot communicate with web clients. Setting the port number to **0** effectively disables the HTTP Server module.

ATTENTION!

The Web IP Port of the HTTP Server must be **forwarded by the router** to allow access from external clients.

Web Server Name

This field is used by the database to create fully qualified study tokens, which are provided to patients for retrieving their studies. If the server does not have a valid DNS name, or if the **Web Server Name** field is left blank, the database interface of the server automatically uses the **Internet IP Address (True IP)** to create a valid fully qualified web address, for example:

[http:// 62.103.211.57:3000/token?4m5RHRCmSrWYcnmk8ICTHvg0qy2KHA7K5ZldW7ZDNsEKsZ68NcJU4b83UtnYJfwx](http://62.103.211.57:3000/token?4m5RHRCmSrWYcnmk8ICTHvg0qy2KHA7K5ZldW7ZDNsEKsZ68NcJU4b83UtnYJfwx)

ATTENTION!

Leave the **Web Server Name** field blank if the server does not have a valid DNS name.

Enable SSL Encryption


Because the server stores sensitive and confidential data, it is **strongly recommended** to use the HTTP Server module with **SSL encryption enabled**, and, if possible, to operate it only within **VPN connections**.

To obtain an SSL certificate, the user must have a valid **ICANN-registered domain name** (for example, *santesoft.com*).

When the **Enable SSL Encryption** option is selected, the program prompts the user to provide two files:

- the **certificate file** (.crt)
- the **private key file** (.key)

These files are issued by the SSL certificate provider chosen by the user.



ATTENTION!

When SSL encryption is enabled, the server can be accessed **only through its domain name** (for example, <https://mycompany.com:3000>) and **not** through its IP address (for example, <https://62.103.211.57:3000> or <https://192.168.1.131:3000>).

Study Accessibility

This combo box allows the user to define the period of time during which a patient or a remote physician can access studies using the token mechanism. The time period is calculated starting from the **study date**. By adjusting this setting, the user can limit the duration of external access and thereby reduce the **server workload** generated by remote users.

How to use the server


The user must use a web browser like **Edge**, **Chrome** or **Mozilla** to communicate with the web server.


If the SSL Encryption is enabled

The user must put in the address bar of the browser (Edge, Chrome or Mozilla) the server's ICAN name followed by a colon and the selected port. The protocol must be the encrypted `https://` protocol.

e.g.

`https://santesoft.gr:3000`



Please note the padlock symbol  before the name of the server. That means that the encrypted protocol is used.

The local IP Address (e.g. 192.168.1.131) or the localhost IP Address 127.0.0.0 do not work because the encrypted `https://` protocol requires a valid server name.

`https://192.168.1.131:3000` (does not work)

`https://127.0.0.1:3000` (does not work)

Furthermore, the simple protocol `http://` (without "s") does not work

`http://santesoft.gr:3000` (does not work)

`http://192.168.1.131:3000` (does not work)

`http://127.0.0.1:3000` (does not work)

If the SSL Encryption is not enabled

The user must enter in the browser's address bar (Edge, Chrome, or Mozilla Firefox) the server's **ICANN-registered domain name**, the **local IP address** (e.g., 192.168.1.131), or the **localhost address** (127.0.0.1), followed by a colon and the selected port number.


The protocol must be the standard `http://` protocol.

e.g.

`http://santesoft.gr:3000`

`http://192.168.1.131:3000`

`http://127.0.0.1:3000`



Furthermore, the encrypted protocol `https://` (with "s") does not work

`https://santesoft.gr:3000` (does not work)

`https://192.168.1.131:3000` (does not work)

`https://127.0.0.1:3000` (does not work)

ATTENTION!

It is strongly discouraged to use the server without SSL encryption. In such cases, it should be operated exclusively within VPN connections.

Edit button

This button is enabled when the user has selected a row of the HTTP Server User list and lets the administrator to add more useful information about the web user like the real name of the web user and his/her contact information.

“Is Active User” option

With this option the administrator can permit or prohibit the login of the user into web server, temporarily or permanently.

“Can upload files” option











If the **“Can upload files”** option is enabled, the query web page of the server displays a section for file uploading. Otherwise, this section does not appear. The files to be uploaded can be DICOM files or zip files that contain DICOM files that will be inserted into PACS server, and text files that contain HL7 messages that will be inserted into Worklist server.

Upload File (a DICOM file, a Zip file that contains DICOM files or a Worklist HL7 text message)	
Step 1: Select a file	Step 2: Upload the file
<input type="button" value="Choose File"/> No file chosen	<input type="button" value="Upload File"/>

“User can burn CD/DVD disks” combo box

It defines if the remote user can burn CD/DVD disks remotely on an Epson or Primera producer.

PATIENT1 FIRST A (67-98-87) , M, 15/03/1975

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
01 (U)	911832595	23/11/1998 14:49:54		US	1	Press the icon to select a PDF file  <input type="button" value="Upload File"/>	         <input type="button" value="Create Study CD/DVD"/>

Change User Name/Password button

By pressing this button, a dialog box appears and lets the user change the **User Name** and the **Password** of a web user. The password must be at least 8 characters and it is **strongly** recommended the rules of selecting a strong password to be followed. That means the password must be a mix of lowercase and uppercase numbers and symbols. This button is enabled when the user has selected a row from the HTTP Server User list.

Sante PACS Server
✕

Change User Name / Password

User Name


Password

Confirm Password

Delete button


This button deletes a web user from the list and it is enabled when the user has selected a row of the HTTP Server User list.

To get a detailed report about every connection in all the three above lists, the user must double-click on the desired row. A dialog box appears and displays all the activity of the web user.




Remote connection with an internet browser

The initial page of the server informs the user that must have an account and must login to connect with the server. Also, it provides links for the two free applications, Sante Web Downloader and Sante DICOM Web Viewer.



Sante Web Downloader lets the user review the downloading procedure. In the end of receiving Sante Web Downloader saves the received studies in a Zip file. The user can configure Sante Web Downloader to open the saved zip file with a DICOM Viewer of his/her choice like Sante DICOM Viewer Pro or Sante DICOM Editor. The DICOM Viewer must support opening zip files from command line. For more information, please read the **User Guide** of Sante Web Downloader. This application is useful when a user wishes to use a DICOM Viewer other than Sante DICOM Web Viewer.

Sante DICOM Web Viewer lets the user review the downloading procedure as Sante Web Downloader does, but also it is a full featured, completely free 2D/3D DICOM Viewer, with all the necessary tools like measurements and study/series comparison and allows the user to work smoothly in a multi-modality, multi-manufacturer environment. For more information, please read the **User Guide** of the application.



Sante Web Downloader



Sante DICOM Web Viewer

The next page of the server asks the remote user for the username and the password. The administrator of the system must create an account for every user of the web server.

The screenshot shows a web browser window displaying the login page of the Santesoft PACS Server. The browser's address bar shows the URL: 192.168.1.6:3000/loginpage?system=santesoft+Ltd. The page features the Santesoft logo and 'Sante PACS Server' branding. Below the header, it indicates 'Santesoft LTD - R&D Department' and 'Client IP: 192.168.1.6'. The main content area contains a 'Login' form with two input fields: 'Username' and 'Password', followed by a blue 'Login' button. The footer of the page contains the copyright notice: 'Copyright © 2022 Santesoft LTD. all rights reserved'.

The next page of the server allows the user to query the PACS server remotely and find the studies of the desired patient and/or the studies of the desired date. The search criteria can be the patient's name, the patient ID, the patient sex, the patient birthdate, the study date, the accession number, the institution, the referring physician and the modality.


From this page the user can also upload files that can be DICOM files and zip files that contain DICOM files that will be inserted into PACS server, and text files that contain HL7 messages that will be inserted into Worklist server. The user can upload files if the administrator of the system has given this right to the particular remote user. Otherwise, the section of the page that allows uploading of files does not appear.

The screenshot displays the query and upload interface of the Santesoft PACS Server. The browser address bar shows the URL: 192.168.1.6:3000/query?session_id=Nq6HjQRKScDs. The page header includes the Santesoft logo and 'Sante PACS Server'. Below the header, it shows 'Santesoft LTD - R&D Department' and 'Client IP: 192.168.1.6'. The main content area is divided into two sections: 'Query' and 'Upload File'. The 'Query' section contains a search form with the following fields: 'Search' (Today), 'Patient Name', 'Patient ID', 'Patient Sex' (All), 'Date of birth' (mm/dd/yyyy), 'Study date from' (mm/dd/yyyy), 'Study date to' (mm/dd/yyyy), 'Accession number', 'Institution', 'Referring Physician', and 'Modality' (ALL). A blue 'Submit' button is located below the query fields. The 'Upload File' section is titled 'Upload File' and includes the description '(a DICOM file, a Zip file that contains DICOM files or a Worklist HL7 text message)'. It has two steps: 'Step 1: Select a file' with a 'Choose File' button (showing 'No file chosen') and 'Step 2: Upload the file' with an 'Upload File' button. The footer of the page contains the copyright notice: 'Copyright © 2022 Santesoft LTD. all rights reserved'.

To search within a range of dates, the user must select the dates in the “Study date from” and “Study date to” and in the “Search” combo box must select the item “With Study Date”.

Query		
	<div style="border: 1px solid red; padding: 2px;"> Search With Study Date ▼ </div>	
Patient Name <input type="text"/>	Patient ID <input type="text"/>	Patient Sex All ▼
Date of birth mm / dd / yyyy <input type="text"/>	<div style="border: 1px solid red; padding: 2px;"> Study date from 01 / 20 / 2022 <input type="text"/> </div>	<div style="border: 1px solid red; padding: 2px;"> Study date to 01 / 31 / 2022 <input type="text"/> </div>
Accession number <input type="text"/>	Institution <input type="text"/>	Referring Physician <input type="text"/>
Modality ALL ▼		
<div style="background-color: #0056b3; color: white; padding: 5px 15px; border-radius: 3px; display: inline-block;">Submit</div>		

After a query the system displays the patient results if any.



There are four buttons for patient result manipulation:














With this button the user can download all the studies of the patient in a zip file.

With this button the user can download all the studies of the patient with Sante Web Downloader.

With this button the user can download all the studies of the patient with Sante DICOM Web Viewer.

With this button the user can review all the studies of the patient.

If the user selects to review all the studies of a patient with the button  a new page appears with the study results.

PATIENT1 FIRST A (67-98-87) , M, 15/03/1975							
#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
01 (U)	911832595	23/11/1998 14:49:54		US	1	Press the icon to select a PDF file  Upload File	        

[Go Back](#)

Press the icon to select a PDF file



Upload File

The remote user can attach a PDF report file to a study by using the **“Upload File”** button of the study’s row. Before that, a PDF file must be selected by pressing the **“Select PDF File”** icon.

There are other five buttons for study result manipulation:



With this button the user can download all the series of the study in a zip file.



With this button the user can download all the series of the study with Sante Web Downloader.



With this button the user can review all the series of the study with Sante DICOM Web Viewer.



With this button the user can review all the series of the study with the OHIF (Open Health Imaging Foundation) Viewer. This viewer works without installation in all systems (PCs, mobile phones, tablets, Windows, MacOS, Linux, Android, iOS etc.)



With this button the user can copy to clipboard the token of the study.



With this button the user can burn a CD/DVD disk with the selected study on remote Epson or Primera producer.



With this button the user can change the Read/Unread status of the selected study.



With this button the user can view the study comment if exists. If the study does not have any comment, this button does not appear.





With this button the user can insert, edit or delete a study comment.






With this button the user can review all the series of the study.



If the user selects to review all the series of a study with the button  a new page appears with the series results.



The screenshot shows a web browser window displaying the Santesoft PACS and DICOM software interface. The page title is "Santesoft PACS and DICOM software" and "Sante PACS Server". The user is logged in as "Santesoft LTD - R&D Department" with a client IP of "192.168.1.6". The main content area displays a table of series results for "Series of Anonymous 1042 (91252), F".

#	Parent Study	Modality	Series Date/Time	Series Description	Files	Action
0001	302	XC	07/08/2022 11:17:52	Fluo L	1	  

Below the table is a "Go Back" button. The footer of the page contains the text "Copyright © 2022 Santesoft LTD, all rights reserved".

Now, there are three buttons for series result manipulation:



With this button the user can download all the files of the series in a zip file.



With this button the user can download all the files of the series with Sante Web Downloader.



With this button the user can review all the series of the study with Sante DICOM Web Viewer.

Study retrieving by patients

Study retrieving with use of token

Patients can retrieve their studies from the web interface of the server. For that purpose, it is used a mechanism of tokens. A **token** is a keyword, 64 characters wide, that is unique for every study and it has a form similar to:

4m5RHRCmSrWYcnmk8lCThvq0qy2KHA7K5ZldW7ZDNsEKsZ68NcJU4b83UtnYJfwx

By using a token like that, patients can have access to their studies from the web interface of the server, without having a user account on the server. The 64 characters wide token is way more secure than a user account with a username usually 8-12 characters wide and a password usually 8-12 characters wide and it does not require extra management like the creation of patient accounts. The patient's page on the server is like that:

The screenshot shows a web browser window displaying the Santesoft PACS and DICOM software interface. The browser address bar shows the URL: `77.49.171.16:3000/token?S3WxByfduCpO5A9oQdr8CnBBBnedEiMhSme...`. The page header includes the Santesoft logo and the text "PACS and DICOM software". The main content area displays the following information:

- Patient name: Anonymous 1042
- Study ID: 302
- Study description:
- Study date/time: 07/08/2022 11:10:22
- Modalities: XC

Below the study details, there are three sections with links:

- Press here to view the study on-line with OHIF Viewer: [View the study on-line](#)
- Press here to view the study on-line with Sante DICOM Web Viewer: [View the study on-line](#)
(Sante DICOM Web Viewer must be installed, [download it from here](#))
- Press here to download the study: [Zip file of the study](#)
- Press here to download a free image viewer for the study: [Sante DICOM Viewer Lite](#)

At the bottom, there are instructions:

1. Download and install Sante DICOM Viewer Lite (requires Windows 64-bit)
2. Download the zip file of the study
3. Run Sante DICOM Web Viewer, select the menu command "File->Open Zip File" and open the zip file of the study

The footer of the page contains the text: "Copyright © 2022 Santesoft LTD, all rights reserved".


This page lets the users download their study and a free program to view them.

Obtaining the token of a study

There are two ways for obtaining the token of the study for sending to patient, from the interface of the database and from the web interface.








From the database

Select the study and press the **“Copy Study Token”** button. The token is copied to clipboard.



From the web interface

Studies of Anonymous 1042 (91252), F

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
01	302	07/08/2022 11:10:22		XC	1	Press the icon to select a PDF file  Upload File	      Copy study's token to clipboard

Go Back

If the browser does not allow copying to clipboard (it is a user option), the user must mark and copy the token manually.



Token with fully qualified web address

To access a study from the web, patients need a token with a fully qualified web address, i.e. a token with the web address of the server behind it: <http://mycompany.com:3000/token?4m5RHRCmSrWYcnmk8lCThvg0qy2KHA7K5ZldW7ZDNsEKsZ68NcJU4b83UtnYJfwx>

The token that is provided from the web interface of the server is always fully qualified, because the web interface knows its web address. The token that is provided by the database interface does not know the web address of the server. To solve this issue the user must set the web address of the server by using the menu command **“Network → HTTP Server Setup”**:


If the server does not have a valid DNS name or the **“Web server name”** field has been kept blank, the database interface of the sever uses the **“Internet IP Address (True IP)”** arithmetic value for creating a token with a valid fully qualified web address:

<http://77.49.171.16:3000/token?4m5RHRCmSrWYcnmk8lCThvg0qy2KHA7K5ZldW7ZDNsEKsZ68NcJU4b83UtnYJfwx>


This fully qualified web address with the token must be sent to the patient. The patient has just to double-click it to go to the study’s page.

Sante DICOM Web Viewer


Sante DICOM Web Viewer is a fully functional, feature rich 2D/3D DICOM Viewer. For more information, please read the **User Guide** of the application.




Frame View – 6 Tiles




Series View – 4 Tiles




MPR View




Scout lines – Localizer




Series Comparison



Built-in 3D Viewer



ROI Statistics



Distance and Angle Measurements

The OHIF (Open Health Imaging Foundation) Viewer

The Open Health Imaging Foundation (OHIF) was launched in 2015 through a development partnership between Massachusetts General Hospital (MGH) Department of Radiology and Radical Imaging to create an open-source web-based medical imaging platform for the global community. The partnership combines the project management, user experience design, and product testing experience of the MGH group with the software development and system design expertise of Radical Imaging. Funding for OHIF has been primarily through research grants to MGH, direct funding from MGH Department of Radiology and its Precision Imaging Metrics clinical trials informatics project, Radical Imaging through in-kind support, and academic-industry technical collaborations among MGH, Radical, and companies using the OHIF Framework. (From the web site of OHIF)

The MIT License

Copyright © 2018 Open Health Imaging Foundation

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

The interface of the OHIF Viewer

The screenshot displays the OHIF Viewer interface within a web browser. The browser tab is titled "OHIF Viewer" and the address bar shows the URL: `192.168.1.6:3000/viewer/1.2.392.200046.100.3.4.300071.302.20220708111310`. The interface includes a top navigation bar with the "Open Health Imaging Foundation" logo, a "Study List" button, and the text "INVESTIGATIONAL USE ONLY" and "Options". Below this is a toolbar with icons for Series, Stack Scroll, Zoom, Levels, Pan, Length, Annotate, Angle, Reset, CINE, More, and Layout. A "Measurements" button is also present. The main viewing area shows a large circular fluoroscopic image of a vessel. In the bottom left corner, there is a thumbnail of the image labeled "FluoL" with "s: 5" and a small "1" icon. The bottom right corner displays technical details: "Zoom: 21%", "W: 207 L: 104", and "Lossy: NaN : 1". The bottom left corner of the main image area shows "Ser: 5" and "4992 x 3328". The top right corner of the main image area shows the date and time: "Jul 8, 2022 11:10:22".

Customize the web page appearance

Select the menu command "Network → HTTP Server Setup".

A dialog box appears, allowing the user to customize the appearance of the pages of the program's built-in website, so that they match the company's colors, logos, and links.

Sante PACS Server
✕

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Update

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.sanl Update

Institution: Santesoft LTD
Department: R&D Department Update
Web Interface Language: English v

Santesoft
PACS and DICOM software
Powered by Santesoft

Home Page
Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Go Back

Copyright © 2017-2025 Santesoft LTD, all rights reserved

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the logo image, the favicon, the web site's address and the popup alternative text of the logo image

Use the highlighted controls below to customize the logo image, the favicon, the web site's address and the popup alternative text of the logo image.

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...


Site's Favicon (16x16 pixels) Select...

Company's Web Site Address Update

Company's Logo Alternative Text Update







Santesoft Software

Institution: Santesoft LTD | Department: R&D Department | Web Interface Language: English









Santesoft LTD - R&D Department Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file  Upload File	    

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file  Upload File	    

[Go Back](#)

Footer's Text (HTML Code) Update

Copyright © 2017-2025 https://www.santesoft.com Copyright © 2017-2025 Santesoft LTD. all rights reserved

OK Restore Previous Values Restore Program's Defaults Cancel

Customize the header's text, text and background color

Use the highlighted controls below to customize the header's text, text and background color.

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Update

Institution: Santesoft LTD Update

Department: R&D Department Update

Web Interface Language: English v

Santesoft

PACS and DICOM software

Powered by Santesoft

Santesoft LTD - R&D Department Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

[Go Back](#)

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.sani Update

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the quote's text and background color

Use the highlighted controls below to customize the quote's text and background color.

Sante PACS Server
✕

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Update

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.san Update

Institution

Department

Update

Web Interface Language: English ▼

PACS and DICOM software
Powered by Santesoft

Home Page
Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Go Back

Copyright © 2017-2025 Santesoft LTD, all rights reserved

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the unread study's text and background color, and the subtitle's text and background color

Use the highlighted controls below to customize the unread study's text and background color, and the subtitle's text and background color.

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Header's Text Update

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.sanl Update

Institution
Santesoft LTD

Department
R&D Department Update

Web Interface Language:
English v

Santesoft
PACS and DICOM software
Powered by Santesoft

Home Page
Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Go Back

Copyright © 2017-2025 Santesoft LTD, all rights reserved

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the read study's text and background color, and the subtitle's text and background color

Use the highlighted controls below to customize the read study's text and background color, and the subtitle's text and background color.

Sante PACS Server ✕

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Header's Text Update

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.san Update

Institution: Santesoft LTD
Department: R&D Department Update
Web Interface Language: English v

PACS and DICOM software
Powered by Santesoft

Home Page Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Go Back

Copyright © 2017-2025 Santesoft LTD, all rights reserved

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the button's text and background color, and the hover text and background color

Use the highlighted controls below to customize the button's text and background color, and the hover text and background color.

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address
https://www.santesoft.com/ Update

Company's Logo Alternative Text
Santesoft Software Update

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Header's Text Update

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code)
Copyright © 2017-2025 <a href="https://www.sani Update

Institution: Santesoft LTD
Department: R&D Department Update
Web Interface Language: English v

Santesoft
PACS and DICOM software
Powered by Santesoft

Home Page Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	Press the icon to select a PDF file Upload File	

Go Back

OK
Restore Previous Values
Restore Program's Defaults
Cancel

Customize the footer's text, text and background color and the company's link

Use the highlighted controls below to customize the footer's text, text and background color and the company's link. The text accepts HTML code.

Web Page Appearance

Company's Logo Image (225x64 pixels) Select...

Site's Favicon (16x16 pixels) Select...

Company's Web Site Address Update

https://www.santesoft.com/

Company's Logo Alternative Text Update

Santesoft Software

#399df5 Header's Text Color Select...

#193fa3 Header's Background Color Select...

Powered by Santesoft Update

Header's Text

#000000 Quote's Text Color Select...

#eef6f9 Quote's Background Color Select...

#000000 Unread Study's Text Color Select...

#98d1ff Unread Study's Background Color Select...

#000000 Read Study's Text Color Select...

#daa8f9 Read Study's Background Color Select...

#000000 Subtitle's Text Color Select...

#cae7a6 Subtitle's Background Color Select...

#ffffff Button's Text Color Select...

#193fa3 Button's Background Color Select...

#ffffff Button's Hover Text Color Select...

#399df5 Button's Hover Background Color Select...

#43bcf1 Footer's Text Color Select...

#193fa3 Footer's Background Color Select...

Footer's Text (HTML Code) Update

Copyright © 2017-2025 <a href="https://www.san

Institution
Department
Web Interface Language:

PACS and DICOM software
Powered by Santesoft

Home Page
Log out
Santesoft LTD - R&D Department
Client IP: 127.0.0.1

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
05 (U)	5685	02/08/2013 19:03:25	THORAX	CT/SR	9	<div style="font-size: 8px;">Press the icon to select a PDF file</div> <div style="font-size: 8px; background-color: #336699; color: white; padding: 2px; text-align: left;">Upload File</div>	

Anonymized (5685), M, 13/12/1999

#	Study ID	Study Date/Time	Study Description	Modalities	Series	Attach PDF To Study	Actions
06 (R)		02/08/2013 19:03:25	THORAX	CT/SR	9	<div style="font-size: 8px;">Press the icon to select a PDF file</div> <div style="font-size: 8px; background-color: #336699; color: white; padding: 2px; text-align: left;">Upload File</div>	

Copyright © 2017-2025 Santesoft LTD, all rights reserved